
**ENGLISH**

**PUSH SPORTS ANKLE BRACE 8**  
Push® is a registered trademark of Net International BV

01/2022  
47204 / 3  
Made in the EU  
Designed in the Netherlands

**PUSH SPORTS**  
**ANKLE BRACE 8**  
Useful information for you!

**WASHING INSTRUCTIONS/MAINTENANCE**  
To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

**DEUTSCH**

**PUSH SPORTS KNÖCKELBANDAGE 8**  
Die Push Sports Ankle Brace 8 provides stability with an individually adjustable elastic strap that wraps around the ankle in the shape of a figure '8'. This substantial bandage is easy to fit and put on. The thin, light material around the foot wicks away moisture, feels comfortable and dries quickly. The brace is easy to use and stays in place well. The brace can be worn without shoes.

**WHEN TO USE THE PUSH SPORTS ANKLE BRACE 8?**  
If you want additional ankle support because you feel insecure or unstable (perhaps due to a previous ligament sprain), the Push Sports Ankle Brace 8 is a good choice. The elastic construction of the brace provides a comfortable level of compression and creates an increased sense of stability.

**INDICATIONS**  

- To prevent recurring ankle injury
- Mild distortion
- Capsule irritation
- Mild instability

**CONTRAINDICATIONS**  

- There are no known disorders where the Push Sports Ankle Brace 8 should not be worn.

**CAUTIONS**  


- Avoid fitting the brace too tightly and be aware of slipping if the brace is worn without a shoe on a slippery surface. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, it is recommended that you fit the brace before purchasing to check its size. When buying online, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

**SPECIAL FEATURES**  
The elastic strap that wraps around the ankle in the shape of a figure '8' stabilises the ankle and provides substantial and comfortable compression. The special TPU dots on the inside of the brace run across the front part of the foot and the heel to prevent the brace from sliding, increasing the effectiveness of the elastic and contact with the skin.

**TAKING MEASUREMENTS**  
The Push Sports Ankle Brace 8 is available in three sizes and in right and left versions. Select the appropriate brace size by measuring the entire circumference of the instep.

**FITTING**  

- Put on the brace. The heel section of the brace needs to fit snugly around the heel.
- Pull your foot towards you, such that the foot is positioned at a 90 degree angle. Then tighten and fasten the top portion of the brace.
- Take the broad elastic strap and wrap it around the back of the lower leg to the front, diagonally across the instep (Fig. ③); then pass it underneath the foot and along the outside of the foot towards the inside of the ankle (Fig. ④). You effectively create a figure '8'. Fasten the strap with the hook and loop material (Fig. ⑤).

[pushsports.eu/ankle-brace-8/applying](https://pushsports.eu/ankle-brace-8/applying)

For details of the distributor in your country, please refer to [www.push.eu/distributors](http://www.push.eu/distributors)

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**De Push Sports Ankle Brace 8 is easy to fit and put on. The thin, light material around the foot wicks away moisture, feels comfortable and dries quickly. The brace is easy to use and stays in place well. The brace can be worn without shoes.**

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