PUSH SPORTS KNEE BRACE

The Push Sports Knee Brace is a very comfortable sports brace that provides knee support with leaf spring hinges along the sides of the knee joint. It prevents lateral movement of the knee joint, without restricting flexion and extension. The knee brace features a U-shaped pad around the patella, which exerts pressure on the patella tendon and stabilizes the patella. The accurate fit and the anti-slip lining prevent the brace from slipping. The thin, light material increases comfort by wicking perspiration away from the skin. Despite its strong support, the Push Sports Knee Brace has a remarkably low profile design, allowing it to be worn under regular trousers. As the leaf spring hinges are never covered, the brace can also be used safely in contact sports. If you have any questions, please consult your retailer or send an e-mail to info@push.eu.

WHEN TO USE THE PUSH SPORTS KNEE BRACE?

Do you need support and comfortable compression around the knee joint during sports activities? Then the Push Sports Knee Brace is a good choice.

The brace is used in the case of residual complaints after ligament injuries or a meniscus surgery. The brace is also effective when the knee is under heavy strain, in sports that require repeated rotation or turning.

INDICATIONS

- General knee instability
- Residual instability following a collateral ligament injury
- Arthritis with mild instability
- Meniscus injury
- Patellofemoral pain syndrome

CONTRAINDICATIONS

- There are no known disorders where the Push Sports Knee Brace cannot or should not be worn.

CAUTIONS

- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports massage therapist.
- The Push Sports Knee Brace provides lateral support, suitable for minor to moderate (residual) instability.
- In some cases, the brace may touch the other leg. Athletes who experience discomfort are advised to wear trousers that cover the location concerned.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, it is recommended that you fit the brace before purchasing to check its size. When buying online, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.

- Check your Push Sports brace each time before use for signs of wear or aging of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

SPECIAL FEATURES

The knee brace provides effective knee support using lightweight, yet very strong bilateral leaf spring hinges, which mimic the natural bending motion of the knee. The use of this type of (non-axial) hinge provides excellent lateral support of the knee joint without causing any restriction during sports activities. The U-shaped pad in the brace surrounds the kneecap and stabilizes the patella.

TAKING MEASUREMENTS

The Push Sports Knee Brace is available in four sizes, universally applicable for the left and right sides. Select the appropriate brace size by measuring the circumference of the lower leg just below the knee. It is recommended to always try out the brace before use, to ensure that it is the right size.

FITTING

Open the two hook and loop fasteners.

- Put the brace over the knee until the pad is positioned snugly underneath the kneecap. For the proper degree of support, the pad must fit closely around the kneecap. Please make sure the lower edge of the brace does not fold inwards.
- Close the upper hook and loop fastener in such a way that you feel a comfortable degree of pressure around the upper leg.
- Close the lower hook and loop fastener in such a way that you feel a comfortable degree of pressure below the knee.

pushsports.eu/knee-brace/applying

For details of the distributor in your country, please refer to www.push.eu/distributors

WASHING INSTRUCTIONS/MAINTENANCE

[🧼](#) [🚿](#) [🧴](#) [🚰](#) [🚰](#) [🚰](#)

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to prevent the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

IMPORTANT

The knee brace is not meant as a replacement for a knee protector. Athletes who use knee protectors can wear the knee protector over the brace.

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