01/2022

47205 / 4

Made in the EU

Designed in the Netherlands

www.push.eu

The Netherlands

6199 AB Maasricht-Airport

Nea International bv

Europalaan 31

PUSH SPORTS

KNEE BRACE

Useful information for you!

Open the two hook and loop fasteners.

1 Pull the brace over the knee until the pad is positioned snugly underneath the kneecap.

For the proper degree of support, the pad must fit closely around the kneecap.

Please make sure the lower edge of the brace does not fold inwards.

2 Close the upper hook and loop fastener in such a way that you feel a comfortable degree of pressure around the upper leg.

3 Close the lower hook and loop fastener in such a way that you feel a comfortable degree of pressure below the knee.

pushsports.eu/knee-brace/applying

For details of the distributor in your country, please refer to www.push.eu/distributors

WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand.

Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

IMPORTANT

The knee brace is not meant as a replacement for a knee protector. Athletes who use knee protectors can wear the knee protector over the brace.

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- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

SPECIAL FEATURES

The knee brace provides effective knee support using lightweight, yet very strong bilateral leaf spring hinges, which mimic the natural bending motion of the knee.

The use of this type of (non-axial) hinge provides excellent medial/lateral support of the knee joint without causing any restriction during sports activities. The U-shaped pad in the brace surrounds the kneecap and stabilizes the patella.

TAKING MEASUREMENTS

The Push Sports Knee Brace is available in four sizes, universally applicable for the left and right sides. Select the appropriate brace size by measuring the circumference of the lower leg just below the knee. It is recommended to always try out the brace before use, to ensure that it is the right size.

FITTING

Open the two hook and loop fasteners.

1 Pull the brace over the knee until the pad is positioned snugly underneath the kneecap.

For the proper degree of support, the pad must fit closely around the kneecap.

Please make sure the lower edge of the brace does not fold inwards.

2 Close the upper hook and loop fastener in such a way that you feel a comfortable degree of pressure around the upper leg.

3 Close the lower hook and loop fastener in such a way that you feel a comfortable degree of pressure below the knee.

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