
Push® is a registered trademark of Nea International bv

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The Push Sports Thumb Brace stabilises, supports and protects the ulnar collateral ligament (MCP-1), leaving the unaffected thumb joints completely free to move. The straps can be adjusted separately, allowing you to determine the most comfortable position for the thumb in your sport. The system's straps cross and provide lateral support to the thumb joint. The Push Sports Thumb Brace leaves the palm of the hand free to hold or make contact with any piece of sports equipment. It fits underneath a skier's or goalkeeper's glove and provides protective support during sports. If you have any questions, please consult your retailer or send an e-mail to info@push.eu.

WHEN TO USE THE PUSH SPORTS THUMB BRACE?

In the event of damage, irritation or pain in the ligaments on the inside of the thumb. This can be the result of hyperextension during ball-sports, or 'skier's thumb'. The brace is also indicated in the event of pain in the entire joint at the base of the thumb. The Push Sports Thumb Brace creates an optimal recovery environment for the ligaments and protects the joints against new injuries.

INDICATIONS

- (Partial) rupture of the ulnar collateral ligament MCP-1 (skier's thumb)
- Distortion of the MCP-1 joint
- Instability in the MCP-1 joint

CONTRAINDICATIONS

- The Push Sports Thumb Brace is not designed to support the CMC-1 joint. In the event of complaints in this joint, such as CMC-1 osteoarthritis, the Push Sports Thumb Brace is not the most suitable solution.

CAUTIONS

- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.
- When using the brace in situations where a glove is typically worn, it is recommended to wear the brace underneath the glove.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, it is recommended that you fit the brace before purchasing to check its size. When buying online, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

SPECIAL FEATURES

The thumb cap and non-elastic strap are held in position by the wrap-around portion at the base of the brace. The thumb cap comfortably encloses the thumb and acts as an anchoring point for the non-elastic strap. The strap in turn restricts the movements of the thumb joint. The level of restriction is individually adaptable.

TAKING MEASUREMENTS

The Push Sports Thumb Brace is available in 3 sizes, and in right and left versions. To determine the right size, measure around the middle of the hand, beneath the knuckles.

FITTING

- Open both hook and loop fasteners.
- For initial use only**
 - Place the thumb cap over the thumb, ensuring that the MCP joint fits in the bowl of the thumb cap. Then close the thumb loop diagonally across the cap, ensuring a snug fit, but avoiding any pinching.
 - Slide the thumb cap across the thumb, ensuring the slightly arched thumb fits snugly in the thumb cap.
 - Close the elastic wrist strap at the base of the brace comfortably.
 - Align and fasten the section of the elastic base with the single blue line to the section on the loop with the single blue line. Tighten the elastic fabric comfortably.
 - Wrap the non-elastic strap around the back of the thumb and align and fasten the strap with the two blue lines on the loop section.
 - Laste, fasten the yellow tab comfortably on the back of the hand.

pushsports.eu/thumb-brace/applying

For details of the distributor in your country, please refer to www.push.eu/distributors

WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

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